Tetrathionate Broth Base acc. to MULLER-KAUFFMANN

For the selective enrichment of salmonellae from various materials, particularly meat, meat products and other foodstuffs.

General Information

This culture medium complies with the recommendations of the DIN Norm 10160 for the examination of meat and the DIN Norm 10181 for the examination of milk.

Mode of Action

Tetrathionate is produced from thiosulfate by adding iodine to the culture medium. Tetrathionate suppresses the growth of coliform and other enteric bacteria. Salmonella, Proteus and some other species of bacteria can reduce tetrathionate and are not inhibited. Calcium carbonate buffers the sulphuric acid, which is liberated when tetrathionate is reduced. Bile promotes the growth of Salmonella, but largely inhibits the accompanying bacteria. Brilliant green suppresses primarily Gram-positive bacteria.

Typical Composition (g/litre)

Meat extract 0.9; peptone from meat 4.5; yeast extract 1.8; sodium chloride 4.5; calcium carbonate 25.0; sodium thiosulfate 40.7; oc bile 4.75.

Also to be added:

Potassium iodide 5.0; iodine 4.0; brilliant green 0.01.

Preparation

Suspend 82 g/litre, heat briefly to the boil and cool rapidly. A sediment of calcium carbonate appears in the turbid broth at the bottom of the tubes.

■ Do not autoclave.

Prior to use add iodine/potassium iodide solution (20 mg/litre) and a 0.1 % solution of brilliant green (10 ml/litre), dispense into test tubes taking care to suspend any precipitate evenly. Avoid further heating.

Preparation of the iodine/potassium iodide solution: Potassium iodide 5 g; iodine 4 g; distilled water 20 ml.

The medium is turbid and green with a white sediment (calcium carbonate).

Quality control

Test strains	Inoculum	Growth after 24 hours
Escherichia coli ATCC 25922	approx. 99 %	≤ 5 %
Salmonella typhimurium ATCC 14028	approx. 1 %	≥ 95 %

Storage

The ready-to-use broth should be prepared and used the same day.

Experimental Procedure and Evaluation

Directly suspend approximately 10 g of sample material in 100 ml Tetrathionate Broth acc. to MULLER-KAUFFMANN.

Incubation: 18-24 hours at 35-37 °C or 43 °C respectively (BÄNFFER 1971, EDEL and KAMPELMACHER 1969).

The resulting cultures should be subjected to further tests.

Literature

BÄNFFER, J.R.: Comparison of the isolation of Salmonellae from human faeces by enrichment at 37 °C and 43 °C. - Zbl. Bakt. I. Orig., 217; 35-40 (1971).

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Ordering Information

Product	Ordering No.	Pack size
Tetrathionate Broth Base acc. to MULLER-KAUFFMANN	1.10863.0500	500 g
Brilliant green (C.I. 42040)	1.01310.0050	50 g
Iodine resublimed	1.04761.0100	100 g